


**Did you know ?**  
**25000+ Distal Radius Volar of Nebula**  
implanted worldwide since 2013.



**NEB - LOCK**  
**System**

**NEBULA SURGICAL PVT. LTD.**



"Nebula" Narmada Park-3, Opp. King's Height, Vidhya Kunj Society Main Road,  
Amin Marg, Rajkot-360005. Gujarat (INDIA) Ph. : +91 281 258 48 32 / 33,



G-1921/3, Almighty Gate Main Road, G.I.D.C. Metoda-360021. Dist :- Rajkot, Gujarat (INDIA)



[www.nebulasurgical.com](http://www.nebulasurgical.com),



[info@nebulasurgical.com](mailto:info@nebulasurgical.com)  
[export@nebulasurgical.com](mailto:export@nebulasurgical.com)



For over 19 years,

Nebula means experience, consistency & reliability proven with time in the field of orthopaedic industry and customer satisfaction. Our aim is to exclusively focus on orthopaedic Implants and surgical Instruments to enhance the surgical techniques which will improve patients post life surgery & will allow surgeons to maximize possibilities of best recovery for patients & live normal life. Nebula is your dedicated one-stop-shop for your all needs of Implants and Instruments.

Our constant efforts, expertise & experience allow us to serve with more than 2000 products like Bone plates & screws, Locking plates & screws, Interlocking nails, Spinal System, Maxillofacial, Arthroscopy implants & Instruments, External Fixator and much more...

‘Supplying to the world’s leading



orthopaedic hospitals’

The ‘NEBULA’ expansion project facility is spread over 6000 Sq. Meters of area with hi-tech multi floor enclosed building with fully automatic A.H.U. (Air Handling Unit) as per norms & guidelines of US FDA (510K) to serve the global requirements. The entire production facility is enclosed in dust-free & highly sterile environment with fully automated latest CNC, VMC, Laser cutting machines of hybrid generation sourced from Japan, USA, Germany, Korea, Sweden. This is most modern infrastructure with state of art manufacturing facilities of this era.





## Straight Locking Plate



### SMALL LOCKING PLATE 2.7 MM (FOR RADIUS / ULNA)

Cat	Length
L105.27	04 HOLE To 10 HOLE



#### Adopt Screw

Locking screw	Cortical screw
2.7 mm	2.7 mm

### SMALL LOCKING PLATE 3.5 MM (FOR RADIUS / ULNA)

Cat	Length
L105	04 HOLE To 14 HOLE
TL105	04 HOLE To 12 HOLE



#### Adopt Screw

Locking screw	Cortical screw
3.5 mm	3.5 mm

### NARROW LOCKING PLATE 4.9 MM (FOR HUMERUS / TIBIA)

Cat	Length
L106	04 HOLE To 16 HOLE
TL106	04 HOLE To 16 HOLE



#### Adopt Screw

Locking screw	Cortical screw
4.9 mm	4.5 mm

### BROAD LOCKING PLATE 4.9 MM (FOR FEMUR)

Cat	Length
L107	05 HOLE To 16 HOLE
TL107	05 HOLE To 16 HOLE



#### Adopt Screw

Locking screw	Cortical screw
4.9 mm	4.5 mm

## Straight Locking Plate



### RECONSTRUCTION LOCKING PLATE 3.5 MM (FOR RADIUS / ULNA)

Cat	Length
L108	04 HOLE To 16 HOLE
TL108	04 HOLE To 12 HOLE



#### Adopt Screw

Locking screw	Cortical screw
3.5 mm	3.5 mm

### RECONSTRUCTION LOCKING PLATE 4.9 MM

Cat	Length
L117	04 HOLE To 16 HOLE



#### Adopt Screw

Locking screw	Cortical screw
4.9 mm	4.5 mm

### 1/3 TUBULAR LOCKING PLATE 3.5 MM (FOR RADIUS / ULNA)

Cat	Length
L120	04 HOLE To 12 HOLE
TL120	04 HOLE To 10 HOLE



#### Adopt Screw

Locking screw	Cortical screw
3.5 mm	3.5 mm

### RCP LOCKING PLATE (ROUND HOLE) 3.5 MM

Cat	Length
L145	04 HOLE To 08 HOLE



#### Adopt Screw

Locking screw
3.5 mm



PROXIMAL HUMERUS LOCKING PLATE

Cat	Length
L113	03 HOLE TO 08 HOLE
TL113	03 HOLE TO 08 HOLE

Adopt Screw

Locking screw	Cortical screw
3.5 mm	3.5 mm



PHILOS LOCKING PLATE

Cat	Length
L123	03 HOLE To 12 HOLE
TL123	03 HOLE To 12 HOLE

Adopt Screw

Locking screw	Cortical screw
3.5 mm	3.5 mm



DISTAL HUMERUS DORSOLATERAL WITH SUPPORT LOCKING PLATE (L/R)

Cat	Length
L129.L/R	03,05,07,09,14 HOLE
TL129.L/R	03,05,07,09,14 HOLE

Adopt Screw

Locking screw	Cortical screw
2.7 & 3.5 mm	3.5 mm



DISTAL HUMERUS DORSOLATERAL WITHOUT SUPPORT LOCKING PLATE (L/R)

Cat	Length
L130.L/R	03,05,07,09,14 HOLE
TL130.L/R	03,05,07,09,14 HOLE

Adopt Screw

Locking screw	Cortical screw
2.7 & 3.5 mm	3.5 mm



DISTAL HUMERUS MEDIAL LOCKING PLATE (L/R)

Cat	Length
L128.L/R	03,05,07,09,14 HOLE
TL128.L/R	03,05,07,09,14 HOLE

Adopt Screw

Locking screw	Cortical screw
2.7 & 3.5 mm	3.5 mm



DISTAL HUMERUS MEDIAL METAPHYSEAL LOCKING PLATE

Cat	Length
L127	07,09,11,13,15 HOLE
TL127	07,09,11,13,15 HOLE

Adopt Screw

Locking screw	Cortical screw
3.5 mm	3.5 mm



EXTRA ARTICULAR DISTAL HUMERUS LOCKING PLATE (L/R)

Cat	Length
L143.L/R	04,06,08,10,12,14 HOLE
TL143.L/R	04,06,08,10,12,14 HOLE

Adopt Screw

Locking screw	Cortical screw
3.5 mm	3.5 mm





## Distal Radius Locking Plate

### Volar Side



#### SMALL T LOCKING PLATE

Cat	Length
L114	03 HOLE To 08 HOLE
TL114	03 HOLE To 08 HOLE

#### Adopt Screw

Locking screw	Cortical screw
3.5 mm	3.5 mm



#### OBLIQUE T LOCKING PLATE

Cat	Length
L109	03 HOLE To 08 HOLE
TL109	03 HOLE To 08 HOLE

#### Adopt Screw

Locking screw	Cortical screw
3.5 mm	3.5 mm



#### DISTAL RADIUS EXTRA ARTICULAR LCP ( HEAD 5 H ) (L/R)

Cat	Length
L158.L/R	03,05 HOLE
TL158.L/R	03,05 HOLE




#### VA DISTAL RADIUS EXTRA ARTICULAR LCP ( HEAD 5 H ) (L/R)

Cat	Length
L158.V.L/R	03, 05 HOLE

#### Adopt Screw

Locking screw	Cortical screw
2.7 mm	2.7 mm



## Distal Radius Locking Plate

### Volar Side




#### DISTAL RADIUS VOLAR LOCKING PLATE (L/R)

Cat	Length
L131.L/R	02 HOLE To 06 HOLE
TL131.L/R	02 HOLE To 06 HOLE

#### Adopt Screw

Locking screw	Cortical screw
2.7 & 3.5 mm	3.5 mm



#### DISTAL RADIUS OBLIQUE ARTICULAR LCP ( HEAD 5 H ) (L/R)

Cat	Length
L159.L/R	03, 05 HOLE
TL159.L/R	03, 05 HOLE

#### VA DISTAL RADIUS OBLIQUE ARTICULAR LCP ( HEAD 5 H ) (L/R)

Cat	Length
L159.V.L/R	03, 05 HOLE

#### Adopt Screw

Locking screw	Cortical screw
2.7 mm	2.7 mm



#### DISTAL RADIUS BUTTRESS LOCKING PLATE (HEAD 5 H) (L/R)

Cat	Length
L160.L/R	03 HOLE
TL160.L/R	03 HOLE




#### VA DISTAL RADIUS BUTTRESS LOCKING PLATE (HEAD 5 H) (L/R)

Cat	Length
L160.V.L/R	03 HOLE

#### Adopt Screw

Locking screw	Cortical screw
2.7 mm	2.7 mm



## Radius Locking Plate

Dorsal Side



### T DISTAL RADIUS DORSAL LOCKING PLATE

Cat	Length
L141	03, 04 HOLE
TL141	03, 04 HOLE

#### Adopt Screw

Locking screw	Cortical screw
2.7 mm	2.7 mm



### L DISTAL RADIUS DORSAL LOCKING PLATE (HEAD 2 H) (L/R)

Cat	Length
L138.L/R	03, 04 HOLE
TL138.L/R	03, 04 HOLE

#### Adopt Screw

Locking screw	Cortical screw
2.7 mm	2.7 mm



### L DISTAL RADIUS DORSAL LOCKING PLATE (HEAD 3 H) (L/R)

Cat	Length
L139.L/R	03, 04 HOLE
TL139.L/R	03, 04 HOLE

#### Adopt Screw

Locking screw	Cortical screw
2.7 mm	2.7 mm



### OBLIQUE DISTAL RADIUS DORSAL LOCKING PLATE (L/R)

Cat	Length
L140.L/R	03, 04 HOLE
TL140.L/R	03, 04 HOLE

#### Adopt Screw

Locking screw	Cortical screw
2.7 mm	2.7 mm



## Radius Locking Plate

Dorsal Side



### STRAIGHT DISTAL RADIUS DORSAL LOCKING PLATE

Cat	Length
L137	05, 06 HOLE
TL137	05, 06 HOLE

#### Adopt Screw

Locking screw	Cortical screw
2.7 mm	2.7 mm



### U DISTAL RADIUS DORSAL LOCKING PLATE (L/R)

Cat	Length
L136.L/R	17 HOLE TOTAL
TL136.L/R	17 HOLE TOTAL

#### Adopt Screw

Locking screw	
2.7 mm	





## Femur Locking Plate

### Proximal Part



#### PROXIMAL FEMORAL LOCKING PLATE (L/R)

Cat	Length
L125.L/R	03 To 11 HOLE

#### Adopt Screw

Locking C.C. screw	C.C. screw	Locking Cannulated screw	Cortical screw	Locking screw
7.3 mm	7.3 mm	4.9 mm	4.5 mm	4.9 mm



#### D.H.S. LOCKING PLATE

Cat	Length
L118.30/35	04, 05, 06, 08, 10, 12 HOLE
TL118.30/35	04, 05, 06, 08, 10, 12 HOLE

#### Adopt Screw

D.H.S. Lag screw	Locking screw	Cortical screw
12.5 mm	4.9 mm	4.5 mm



#### D.H.S. TROCHANTERIC SUPPORT LOCKING PLATE

Cat	Angle	Length
L118.T.130/135	130°, 135°	04, 05, 06 HOLE

#### Adopt Screw

D.H.S. Lag screw	Cortical screw	Locking screw
12.5 mm	4.5 mm	4.9 mm



## Femur Locking Plate

### Distal Part



#### DISTAL FEMORAL LOCKING PLATE (L/R)

Cat	Length
L112.L/R	04 HOLE To 14 HOLE
TL112.L/R	05, 07, 09, 11, 13 HOLE

#### Adopt Screw

Locking Cancellous screw	Locking screw	Cortical screw
6.5 mm	4.9 mm	4.5 mm



#### 95° D.C.S. LOCKING PLATE

Cat	Length
L119	06, 08, 10, 12 HOLE
TL119	06, 08, 10, 12 HOLE

#### Adopt Screw

D.H.S. Lag screw	Locking screw	Cortical screw	C.C. & Cancellous screw
12.5 mm	4.9 mm	4.5 mm	6.5 mm



#### ANTERO MEDIAL DISTAL FEMUR LCP (L/R)

Cat	Length
L165.L/R	02, 04, 06, 08 HOLE
TL165.L/R	02, 04, 06, 08 HOLE

#### Adopt Screw

Locking screw	Cortical screw
4.9 mm	4.5 mm



## Tibia Locking Plate

### Proximal Part



#### LATERAL TIBIAL LOCKING PLATE (L/R)

Cat	Length
L116.L/R	03 HOLE To 14 HOLE
TL116.L/R	03 HOLE To 11 HOLE

#### Adopt Screw

Locking Cancellous screw	Locking screw	Cortical screw
6.5 mm	4.9 mm	4.5 mm



#### ANATOMICAL PROXIMAL TIBIA LOCKING PLATE (L/R)

Cat	Length
L132.L/R	04, 06, 08, 10, 12, 14 HOLE

#### Adopt Screw

Locking screw	Cortical screw
3.5 & 4.9 mm	4.5 mm



#### ANATOMICAL LOW BEND PROXIMAL TIBIA LCP (L/R)

Cat	Length
L162.L/R	04, 06, 08, 10, 12, 14 HOLE
TL162.L/R	04, 06, 08, 10, 12, 14 HOLE

#### Adopt Screw

Locking screw	Cortical screw
3.5 & 4.9 mm	4.5 mm



#### PROXIMAL TIBIA MEDIAL LCP (L/R)

Cat	Length
L144.L/R	04,06,08,10,12,14,16,18 HOLE
TL144.L/R	04,06,08,10,12,14,16,18 HOLE

#### Adopt Screw

Locking screw	Cortical screw
3.5 mm	3.5 mm



## Tibia Locking Plate

### Proximal Part



#### POSTERIOR MEDIAL PROXIMAL TIBIAL LOCKING PLATE

Cat	Length
L133	02, 04, 06, 08, 10 HOLE
TL133	02, 04, 06, 08, 10 HOLE

#### Adopt Screw

Locking screw	Cortical screw
3.5 mm	3.5 mm



#### L BUTTRESS LOCKING PLATE (L/R)

Cat	Length
L110.L/R	04 HOLE To 12 HOLE

#### Adopt Screw

Locking screw	Cortical screw	C.C. & Cancellous screw
4.9 mm	4.5 mm	6.5 mm



#### T BUTTRESS LOCKING PLATE

Cat	Length
L111	04 HOLE To 12 HOLE
TL111	04 HOLE To 12 HOLE

#### Adopt Screw

Locking screw	Cortical screw	C.C. & Cancellous screw
4.9 mm	4.5 mm	6.5 mm





## Tibia Locking Plate

### Distal Part



#### DISTAL TIBIAL MEDIAL LOCKING PLATE (L/R)

Cat	Length
L121.L/R	04, 06, 08, 10, 12 HOLE
TL121.L/R	04, 06, 08, 10, 12 HOLE

#### Adopt Screw

Locking screw	Cortical screw
3.5 & 4.9 mm	3.5 & 4.5 mm



#### DISTAL TIBIA METAPHYSEAL LOCKING PLATE

Cat	Length
L115	04 HOLE To 12,14,16 HOLE
TL115	04 HOLE To 12 HOLE

#### Adopt Screw

Locking screw	Cortical screw
3.5 & 4.9 mm	3.5 & 4.5 mm




#### ANTEROLATERAL DISTAL TIBIA LOCKING PLATE (L/R)

Cat	Length
L134.L/R	04, 06, 08, 10, 12 HOLE
TL134.L/R	04, 06, 08, 10, 12 HOLE

#### Adopt Screw

Locking screw	Cortical screw
3.5 & 4.9 mm	3.5 & 4.5 mm



## Clavicle Locking Plate



#### SUPERIOR ANTERIOR CLAVICLE LOCKING PLATE (L/R)

Cat	Length
L126.L/R	03 HOLE To 08 HOLE
TL126.L/R	03 HOLE To 08 HOLE

#### Adopt Screw

Locking screw	Cortical screw
2.7 & 3.5 mm	3.5 mm



#### CLAVICLE MEDIAL LOCKING PLATE (L/R)

Cat	Length
L122.L/R	06 HOLE To 08 HOLE
TL122.L/R	06 HOLE To 08 HOLE

#### Adopt Screw

Locking screw	Cortical screw
3.5 mm	3.5 mm



#### CLAVICLE HOOK LOCKING PLATE (L/R)

Cat	Length
L142.L/R	04 HOLE To 08 HOLE
TL142.L/R	04 HOLE To 08 HOLE

#### Adopt Screw

Locking screw	Cortical screw
3.5 mm	3.5 mm



## Fibula

### Distal Part



#### DISTAL FIBULA VA LOCKING PLATE (W.O.S) (L/R)


Cat	Length
L161.W.L/R	04, 06, 08, 10 HOLE
TL161.W.L/R	04, 06, 08, 10 HOLE

#### DISTAL FIBULA VA LOCKING PLATE (L/R)

Cat	Length
L161.L/R	04, 06, 08, 10 HOLE

#### Adopt Screw

Locking screw	Cortical screw
2.7 & 3.5 mm	3.5 mm



## Ulna


### Proximal Part

#### OLECRANON VA LOCKING PLATE

Cat	Length
L135	07,09,11 HOLE (SMALL & LARGE) (L/R)
TL135	07,09,11 HOLE (SMALL & LARGE) (L/R)

#### Adopt Screw

Locking screw	Cortical screw
2.7 & 3.5 mm	3.5 mm



5 - Hole is universal. (No Left or Right)


## Calcaneal

#### CALCANEAL LOCKING PLATE (L/R)

Cat	Length
L124.L/R	SMALL & MEDIUM
TL124.L/R	SMALL & MEDIUM

#### Adopt Screw

Locking screw
3.5 mm



## Locking Screw



#### LOCKING SCREW 2.7 MM

Cat	Length
L152	12 mm To 60 mm (2 mm variation)
TL152	12 mm To 60 mm (2 mm variation)

#### Adopt Instrument

Drill Bit :- 2.0 mm

#### LOCKING SCREW 3.5 MM

Cat	Length
L101	10 mm To 50 mm (2 mm variation) 55 mm To 90 mm (5 mm variation)
TL101	10 mm To 50 mm (2 mm variation) 55 mm To 90 mm (5 mm variation)

#### Adopt Instrument




Drill Bit :- 2.8 mm

#### LOCKING SCREW 4.9 MM

Cat	Length
L102	16 mm To 60 mm (2 mm variation) 65 mm To 120 mm (5 mm variation)
TL102	20 mm To 50 mm (2 mm variation) 55 mm To 100 mm (5 mm variation)

#### Adopt Instrument

Drill Bit :- 4.3 mm





## Locking Cancellous Screw



### LOCKING CANCELLOUS SCREW 3.5 MM (16 THREAD)

Cat	Length
L153	18 mm To 30 mm (2 mm variation) 35 mm To 70 mm (5 mm variation)
TL153	18 mm To 22 mm (2 mm variation) 25 mm To 70 mm (5 mm variation)

#### Adopt Instrument

Drill Bit :- 2.8 mm



### LOCKING CANCELLOUS SCREW 3.5 MM (FULLY THREAD)

Cat	Length
L154	18 mm To 22 mm (2 mm variation) 25 mm To 90 mm (5 mm variation)
TL154	18 mm To 22 mm (2 mm variation) 25 mm To 90 mm (5 mm variation)

#### Adopt Instrument

Drill Bit :- 2.8 mm



### LOCKING CANCELLOUS SCREW 4.9 MM (FULLY THREAD)

Cat	Length
L150	50 mm To 90 mm (5 mm variation)

#### Adopt Instrument

Drill Bit :- 4.3 mm



### LOCKING CANCELLOUS SCREW 6.5 MM (32 THREAD)

Cat	Length
L104	50 mm To 110 mm (5 mm variation)

#### Adopt Instrument

Drill Bit :- 4.5 mm




### LOCKING CANCELLOUS SCREW 6.5 MM (FULLY THREAD)

Cat	Length
L103	35 mm To 110 mm (5 mm variation)
TL103	35 mm To 110 mm (5 mm variation)

#### Adopt Instrument

Drill Bit :- 4.5 mm



## Locking C.C. Screw



### LOCKING C.C. SCREW 7.3 MM (FULLY THREAD)

Cat	Length
L151	50 mm To 151 mm (5 mm variation)

#### Adopt Instrument

Drill Bit :- 5.0 mm C.C.



## Locking Cannulated Screw

### LOCKING CANNULATED SCREW 4.9 MM

Cat	Length
L157	70 mm To 90 mm (5 mm variation)

#### Adopt Instrument

Drill Bit :- 4.3 mm C.C.





1 - Proximal Humerus Locking Plate Graphic Box

Cat No. :- GB018


Length :- 13"  
Width :- 8.5"  
Height :- 3"



2 - Humerus Locking Plate Graphic Box

Cat No. :- GB011

Length :- 14.5"  
Width :- 10"  
Height :- 5.5"



3 - Olecranon Locking Plate Graphic Box

Cat No. :- GB014


Length :- 12"  
Width :- 6"  
Height :- 4"



4 - Distal Radius Locking Plate Graphic Box

Cat No. :- GB012

Length :- 12"  
Width :- 6"  
Height :- 4"









5 - Proximal Tibia Locking Plate Graphic Box

Cat No. :- GB016  
 Length :- 14.5”  
 Width :- 11”  
 Height :- 5.5”




6 - Distal Tibia Locking Plate Graphic Box

Cat No. :- GB017  
 Length :- 14.5”  
 Width :- 11”  
 Height :- 5.5”



7 - Fibula Locking Plate Graphic Box

Cat No. :- GB010  
 Length :- 12”  
 Width :- 6”  
 Height :- 4”



8 - Clavical Locking Plate Graphic Box

Cat No. :- GB007  
 Length :- 12”  
 Width :- 6”  
 Height :- 4”

